

# gBlocks® Gene Fragments

## Protocol: Fragment Amplification



gBlocks Gene Fragments are normalized to 200 ng prior to delivery, which is a suitable quantity for many applications. However, in some cases it may be necessary or desired to amplify your gBlocks Gene Fragments to get additional starting material.

The following information will help you amplify your gBlocks Gene Fragments using the high-fidelity Phusion® DNA Polymerase (www.NEB.com/Phusion). Only use high-fidelity, proofreading enzymes to amplify your gBlocks Gene Fragment to limit the introduction of sequence errors.

### gBlocks Gene Fragments

gBlocks Gene Fragments are chemically synthesized, double-stranded DNA that are compatible with a wide range of existing applications that require double-stranded DNA. gBlocks Gene Fragments are normalized to 200 ng and delivered dried down. Order at [www.idtdna.com/gblocks](http://www.idtdna.com/gblocks).

### Resuspending your gBlocks Gene Fragments

The dried down gBlocks Gene Fragment pellet can become displaced from the bottom of the tube during shipping. Follow these instructions to properly resuspend the pellet for a final concentration of 10 ng/μL.

1. Centrifuge the tube for 3–5 sec at a minimum of 3000 x g to pellet the material to the bottom of the tube.
2. Add 20 μL TE to the tube for a final concentration of 10 ng/μL.
3. Briefly vortex and centrifuge.

### Storing your gBlocks Gene Fragments

gBlocks Gene Fragments can be stored in TE at –20°C for up to 24 months. If gBlocks Gene Fragments will be stored for less than 1 month, they can be resuspended in nuclease-free water instead of TE.

### Designing primers for amplification

For best results, we recommend designing amplification primers for your gBlocks Gene Fragment sequence using the PrimerQuest® tool ([www.idtdna.com/PrimerQuest](http://www.idtdna.com/PrimerQuest)).

### Required materials

- gBlocks Gene Fragments
- Amplification primers
- Phusion® DNA Polymerase ([www.NEB.com/Phusion](http://www.NEB.com/Phusion))
- PCR or gel purification kit (depending on application)

### Phusion® DNA Polymerase amplification reaction

1. Set up the amplification reaction on ice—reaction components for 20 μL and 50 μL reactions are shown.

| gBlocks® Gene Fragments Amplification |                       |                       |
|---------------------------------------|-----------------------|-----------------------|
|                                       | 20 μL Reaction        | 50 μL Reaction        |
| Nuclease-free H <sub>2</sub> O        | Adjust to final 20 μL | Adjust to final 50 μL |
| 5X Phusion HF or GC Buffer            | 4 μL                  | 10 μL                 |
| 10 mM dNTPs                           | 0.4 μL                | 1 μL                  |
| 10 μM Forward Primer                  | 1 μL                  | 2.5 μL                |
| 10 μM Reverse Primer                  | 1 μL                  | 2.5 μL                |
| gBlocks® Gene Fragments               | 0.1–1.0 ng            | 0.1–1.0 ng            |
| Phusion® DNA Polymerase               | 0.2 μL                | 0.5 μL                |
| Total volume                          | 20 μL                 | 50 μL                 |

Table source: [www.NEB.com/Phusion](http://www.NEB.com/Phusion).

2. Gently mix the reaction and spin down in microcentrifuge.
3. Carry out the amplification reaction in a thermocycler with heated lid.

### General cycling conditions

The table shows general guidelines for amplifying a short, 500 bp gBlocks Gene Fragment. Conditions will vary depending on sequence of the gene fragment and amplification primers.

| gBlocks® Gene Fragments Amplification |                    |             |               |
|---------------------------------------|--------------------|-------------|---------------|
| Step                                  | Cycles             | Temperature | Time          |
| Initial denaturation                  | 1                  | 98°C        | 30 seconds    |
| Denaturation                          | 15–25 <sup>1</sup> | 98°C        | 10 seconds    |
| Annealing                             |                    | 45–72°C     | 10–30 seconds |
| Extension                             |                    | 72°C        | 15–30 seconds |
| Final extension                       | 1                  | 72°C        | 5 minutes     |
| Hold                                  | 1                  | 4°C         | ∞             |

<sup>1</sup>Use the lowest number of cycles to attain the desired quantity as this will provide the best sequence fidelity.

Table source: [www.NEB.com/Phusion](http://www.NEB.com/Phusion).

### Purification

Gel purification using high-energy, UV light damages DNA, negatively impacting downstream application performance. We recommend that you purify your gBlocks Gene Fragment amplification product using a non-UV method. Examples of purification methods that do not rely on high-energy, UV light include affinity-column purification, or non-UV gel systems FlashGel™ System (Lonza).